

## Evaluation of cough assist devices by the Swiss Society of Pulmonology (SIG Ventilation & O2)

**Cough Assist Device:** Nippy Clearway 2

**Reviewer / institution:** SSP-summary of evaluation by Georg Volk and team, certified respiratory therapist (USZ)

Device overview	
<b>Company</b>	Breas Medical
<b>Mechanical insufflator exsufflator</b>	NIPPY Clearway 2
<b>CE marking</b>	2797
<b>Product class</b>	<input checked="" type="checkbox"/> Mechanical insufflator exsufflator <input checked="" type="checkbox"/> other: IPPB, NIV (only for intermittent therapeutic ventilation, max. 15 minutes)
Device specifications	
<b>Weight, size</b>	Weight: 3.75 kg (without internal battery) Size/dimension: 285 x 285 x 195 mm
<b>Display</b>	Touchscreen: <input checked="" type="checkbox"/> YES <input type="checkbox"/> NO Size/dimension: 150 x 90 mm Option to dim <input type="checkbox"/> YES <input checked="" type="checkbox"/> NO
<b>Real-time and retrospective monitoring on ventilator display</b>	<p><b>Real-time monitoring</b></p> Curves: <input checked="" type="checkbox"/> pressure-time <input type="checkbox"/> volume-time <input type="checkbox"/> flow-time Scales: <input checked="" type="checkbox"/> pressure <input type="checkbox"/> volume <input checked="" type="checkbox"/> flow  <input checked="" type="checkbox"/> Tidal volume <input checked="" type="checkbox"/> Pressure; specify: cmH2O, hPa <input checked="" type="checkbox"/> CPF <input type="checkbox"/> Te <input type="checkbox"/> Rise time / ramp; specify: <input type="checkbox"/> Other: <input checked="" type="checkbox"/> Supplementary external parameters: <input type="checkbox"/> CO2, if yes <input type="checkbox"/> ETCO2 <input type="checkbox"/> tcpCO2 <input checked="" type="checkbox"/> SpO2 <input type="checkbox"/> other: Adaptable time period for monitoring: <input type="checkbox"/> YES <input checked="" type="checkbox"/> NO if YES, time window: -  Comment: In M-I/E mode, a volume is displayed during inspiration in L. Value difficult to evaluate. NIV mode: Setting corresponds to an APCV mode, only pressure curve can be displayed, pressure and flow additionally on a scale. No display of Vt possible. O2 administration at the expiratory adapter possible.
	<b>Retrospective monitoring (previous usage)</b>

	<input checked="" type="checkbox"/> YES <input type="checkbox"/> NO if YES: <input checked="" type="checkbox"/> Utilisation <input checked="" type="checkbox"/> Profile <input checked="" type="checkbox"/> PCF <input type="checkbox"/> Other: <input type="checkbox"/> Curves: <input type="checkbox"/> Pressure/time <input type="checkbox"/> Volume/time <input type="checkbox"/> Flow/time <input type="checkbox"/> Supplementary external parameters: <input type="checkbox"/> CO2 <input type="checkbox"/> SpO2 <input type="checkbox"/> FiO2 <input type="checkbox"/> other:  Available time periods: 30 + 90 days
<b>Connections</b>	<input checked="" type="checkbox"/> Power adapter, voltage: 100-240V <input checked="" type="checkbox"/> External battery: <input type="checkbox"/> Oxygen inlet port <input type="checkbox"/> Outlet port (for tube) <input checked="" type="checkbox"/> Supplementary measurements/monitoring, e.g. SpO2, CO2, FiO2, p <sub>oes</sub> ): SpO2 <input checked="" type="checkbox"/> Remote control <input type="checkbox"/> Other:  Comment: O2 can only be supplied via expiration valve with IPPB and NIV
<b>Electrical rating and system interface</b>	AC (alternating current): 100-240 VAC, max. 260 W Accu level indicator: <input checked="" type="checkbox"/> YES <input type="checkbox"/> NO
<b>Service</b>	Emergency service (24h/7d): <input checked="" type="checkbox"/> YES <input type="checkbox"/> NO Company: VitalAire Suppliers of disposable material: VitalAire and others
<b>Ventilation modes and settings</b>	
<b>Permitted age groups</b>	Adults <input checked="" type="checkbox"/> YES <input type="checkbox"/> NO Children <input checked="" type="checkbox"/> YES <input type="checkbox"/> NO (not tested in children by SGP)
<b>Device modes</b>	<input checked="" type="checkbox"/> M-I/E <input checked="" type="checkbox"/> IPPB <input checked="" type="checkbox"/> NIV <input type="checkbox"/> Other:
<b>Profiles</b>	Number: 4
<b>Pressure and flow range</b>	M-I/E: Insufflation 3 - 70 cmH2O, Exsufflation -3 - -70 cmH2O IPPB: Max. Insp. P 3 – 60 cmH2O, Flow: 20 – 120 L/min NIV: IPAP 3 – 50 cmH2O Maximum air flow: 120 L/min
<b>Trigger</b>	Inspiratory trigger: <input type="checkbox"/> scale <input checked="" type="checkbox"/> absolute value - Type: <input checked="" type="checkbox"/> flow trigger <input type="checkbox"/> pressure trigger - Range: 0 - 10 - Clarity of the scale: <input type="checkbox"/> arbitrary <input checked="" type="checkbox"/> clear/distinct Comment: No expiratory trigger
<b>Pressure rise</b>	Type: <input checked="" type="checkbox"/> scale <input type="checkbox"/> time Range: 1 - 10 Clarity of the scale: <input type="checkbox"/> arbitrary <input checked="" type="checkbox"/> clear/distinct
<b>Inspiratory and expiratory time</b>	T insp.: NIV 0.5-3.0 s, M-I/E 0.5 – 5.0 s

	T exp: M-I/E 0.5 – 5.0 s Pause: Off / 0.1 – 9.0 s (depends on mode)
<b>Oszillation</b>	Insuff.: frequency 0 – 20 Hz, amplitude 0 – 10 cmH2O Exsuff: frequency 0 – 20 Hz, amplitude 0 – 10 cmH2O
<b>Oxygen</b>	Direct O2 connection on the ventilator: <input type="checkbox"/> YES <input checked="" type="checkbox"/> NO
<b>Alarms</b>	Pressure high/low <input checked="" type="checkbox"/> YES <input type="checkbox"/> NO <input type="checkbox"/> Range: 120% > 2 s Device in battery mode <input type="checkbox"/> YES <input type="checkbox"/> NO Other: No exhalation valve detected, Internal temperature high, Battery level low, Leakage low, Device error, NIV treatment about to end, NIV treatment finished, Power failure, Pressure measurement error, Pressure sensor error, SpO2 sensor disconnected, Treatment too long (to be stored inTreatRepeat), Treatment finished. Option to deactivate all alarms <input type="checkbox"/> YES <input checked="" type="checkbox"/> NO
<b>Safety</b>	Expert area enabling: YES Battery capacity shown <input checked="" type="checkbox"/> YES <input type="checkbox"/> NO If yes: scale: in %
<b>Additions</b>	
<b>Filter</b>	Filter class: Reusable filter (grey), disposable filter (white)
<b>Relevant accessories</b>	Pedestal/Carrier <input type="checkbox"/> YES <input checked="" type="checkbox"/> NO External battery <input checked="" type="checkbox"/> YES <input type="checkbox"/> NO SpO2 sensor <input checked="" type="checkbox"/> YES <input type="checkbox"/> NO Other: Remote control, foot pedal
<b>Software</b>	
<b>Connections / network</b>	<input checked="" type="checkbox"/> USB <input checked="" type="checkbox"/> SD Card, memory size: <input type="checkbox"/> Other:
<b>Software / statistics</b>	Comment: not tested
<b>Handbook</b>	
<b>Manual</b>	Practicability: good Completeness: given
<b>Validation studies</b>	
<b>Validation studies</b>	Not available

<b>Patient report</b>	
<b>Pros / advantages</b>	<u>M-I/E:</u> Test performed on clients with Muscular Dystrophy Type Duchenne both experienced and inexperienced in cough assist use. All reported comfortable use. The therapy was comparable to the experience with other cough assist devices. The operation via the remote control could be carried out independently by patients after brief instruction, despite very limited finger motor skills. Patients found this feature very pleasant.  <u>IPPB:</u> Subjectively, the patient succeeds in achieving complete distribution. The therapy was not perceived as unpleasant with adjusted settings. The patients were motivated to perform the therapy independently after instruction. From a therapeutic point of view, the device proved to be efficient in recruitment and mobilisation of secretions.
<b>Cons / drawbacks</b>	

<b>Summary</b>	
<b>Pros / advantages</b>	<p><u>M-I/E:</u> Easy and clear operation. The application proves to be equivalent overall compared to other devices. In the test of the application via mask and tracheostoma, the Nippy fulfilled the expectations. The patients did not express any undesirable side effects that could not have been remedied by changing the settings. Remote control very practical, can possibly also be operated by the patient himself. TreatRepeat function facilitates the adjustment of the M-I/E</p> <p><u>IPPB:</u> Simple and clear operation. The "plateau function" in particular allows great therapeutic flexibility. The Nippy also proves easy to use for the patient. In the test of use via mask and tracheostoma, the Nippy fulfilled the expectations from a therapeutic point of view.</p>
<b>Cons / drawbacks</b>	<p><u>NIV mode:</u> Setting corresponds to an APCV mode, only pressure curve can be displayed, pressure and flow additionally on one scale. No display of the Vt possible. O2 administration at the expiratory adapter possible. Tube change necessary to switch from M-I/E to NIV. During the test phase, the patient selection did not indicate the use of the additional NIV mode.</p> <p><u>M-I/E:</u> there is no possibility to have each insufflation triggered by the patient. Pat. can only trigger the first breath in basic auto mode. A Vin is measured, which is given in L and is difficult to evaluate. According to the manual, the estimated lung volume during insufflation is given in ml.</p> <p>Despite the set language German in the Compliance Overview and Alarm Log menu, the display is partially in English.</p>
<b>Recommended setting and patient groups</b>	<p>According to manufacturer: <u>M-I/W:</u> muscular dystrophy, amyotrophic lateral sclerosis (ALS) or motor neurone disease (MND), myasthenia gravis, poliomyelitis, other neurological disorders, spinal cord injury. <u>IPPB:</u> Cystic fibrosis, bronchiecatsia, emphysema, post-operative patients (without contraindications).</p>

### **Recommendation of the SIG ventilation & O2 of the Swiss Society of Pulmonology**

Based on the above evaluation, we recommend the availability of this device on the Swiss market.

Comment: IPPB and NIV mode more suitable in clinical than in out-of-hospital settings

23.09.2022

PD Dr. med. Esther I Schwarz



On behalf of the SIG Ventilation and O2